

PILODIST® 107



www.pilodist.de

Tel: ++ 49 (0) 2225 955910

Fax: ++ 49 (0) 2225 9559111

e-mail : info@pilodist.de

PILODIST GmbH
Eichelnkampstraße 2
D-53340 Meckenheim
Sitz: Bonn

Eingetragen AG Bonn HRB 12941
Managing Director Rolf Jeub
USt.-ID.-Nr. DE 813 408 398

Sparkasse KölnBonn
Kto. 33302365
BLZ 37050198
IBAN DE 05 3705 0198 0033 3023 65
Swift-BIC: COLSDE33

PILODIST® 107

PILODIST 107 is an easy to handle solvent recovery system for extremely pure solvents.

The solvent mixtures are side products of extraction's, washing processes, chemical and biochemical reactions, chromatographic operations in clinical range, especially in anatomical and pathological institutes.

The system is equipped with a Concentric-tube-column, i.e. highest separation efficiency and high load ranges.

Characteristics

- Trouble-free operation due to static all glass construction
- All-glass construction
- System is mounted in a safety cabinet with transparent front door.

Application:

All kind of solvents like acetone, DMF, THF, paraffins, Xylenes, alcohol's, acetonitriles, isooctane, hexane, methyl chloride and other HPLC-solvents.

Technical Data

Flask size:	10 l or 20 l
Separation efficiency:	up to 60 theoretical plates
Separation length:	500 mm
Load range:	approx. 150-3000 ml/h
Operating temperature:	max. 300° C
Operating pressure:	atmospheric (optional vacuum)
Material:	borosilikat glass, grease-free
Dimensions (w x h x d):	1.00 x 2.10 x 0.60 m

www.pilodist.de

Tel: ++ 49 (0) 2225 955910

Fax: ++ 49 (0) 2225 9559111

e-mail : info@pilodist.de

PILODIST GmbH
Eichelnkampstraße 2
D-53340 Meckenheim
Sitz: Bonn

Eingetragen AG Bonn HRB 12941
Managing Director Rolf Jeub
USt.-ID.-Nr. DE 813 408 398

Sparkasse KölnBonn
Kto. 33302365
BLZ 37050198
IBAN DE 05 3705 0198 0033 3023 65
Swift-BIC: COLSDE33